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
Oclean Smart Sonic Electric Toothbrush User Manual

Low battery reminder: The toothbrush will not work when battery is lower than 5%, short press on the power button, red light flashes twice and motor vibrates twice to remind that battery is extremely low, recharge is needed in this case.

The red indicator flashes when battery is low. The white indicator is on constantly when battery is full.

Specifications

Product Name	Oclean Sonic Electric Toothbrush
Model	Oclean F1
Max Torsion	140 gf.cm
Waterproof	IPX7
Power	1.8W
Charger	USB charger
Input	5V =0.6A
Battery Life *	30 days
Battery Capacity	800 mAh
Charging Time	2 hours

* Calculated by 2 minutes/time, twice/day.

Facebook Page

www.facebook.com/ocleanofficial

Follow the Oclean Official Facebook Page

Oclean

Remark

- The content of toxic and harmful substance in all homogeneous materials of this component is less than the limited quantity as specified in GB/T 26572-2011;
- × The content of toxic and harmful substance in at least one of the homogeneous materials of this component is over the limited quantity as specified in GB/T 26572-2011;
- * Most materials of this product adopt environment-friendly materials. You can use this product without worry of infringing environment-protection policies.

Name of parts	Toxic and harmful substance or element						Environment Protection Term
	Pb	Hg	Cd	Cr6+	PBB	PBDE	
Brush head	○	○	○	○	○	○	Ten Years
Handle	×	○	○	○	○	○	
Charging base	×	○	○	○	○	○	
Accessory	×	○	○	○	○	○	Five Years
Battery	○	○	○	○	○	○	

Oclean

Warranty Card

- These clauses only apply to Oclean products sold in the country or region where the purchase is made, if there is any inconsistency or omission with the national policy, the national policy shall prevail.
- For performance issues caused by the product itself, one may avail service as per below regulations with valid proof of purchase:
- 2.1 A customer can avail warranty service for free within 24 months from the date of purchase if there is a performance failure caused by the product itself.
- In order to avail replacement and warranty services, one should show the valid proof of purchase.
- Other product related materials (such as Packaging box, Product Manual etc.) are not covered in warranty.

Serial Number:		Name:				
Product model:		Teli:				
		Email:				
		Address:				
		Date of purchase:				
		Invoice number:				
		Problem Description:				