**Go
the
extra smile**

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Important Safety Tips

This device is recommended for use by children aged 6 years and above. For the first use, it is recommended that parents supervise and guide the children for quicker familiarity of the device. When using for the first time, there may be feelings of stimulation or tingling. Please do not worry since it is a normal phenomenon, which will gradually disappear. It is normal if the main body of the device is slightly warm when charging. If the device turns hot, please cut off the power supply immediately and contact the after-sales center. The device can only be used with specified brush heads. Oclean will NOT bear any consequences caused by non-official or authenticated brush heads. The device can only be used for cleaning teeth, gingiva and tongue. Do NOT apply to other usages. The device is for personal cleaning only and is not recommended for medical purpose or multi-user use. The device can only be used indoors. Do NOT place the device under the blazing sun for a long time. It shall be placed in a dry and well-ventilated place with normal temperature. Please cut off the power supply when not charging or before cleaning. Clean the power supply socket at regular intervals to prevent from dusts and dirt accumulation and avoid reducing the charging effect. This device complies with safety standards for electromagnetic devices. For any questions involving pacemakers or other implanted devices, please consult a physician or the provider of the aforementioned device before use. Do not immerse the charging cable in water or any liquid. Do not connect to the power supply when there is liquid in the power plug of the charging cable. With lithium battery inside, the device should be stored away from fire and heat source. Do NOT charge, transport, use and place it in high temperature. Do not knock or drop the device. If the device is damaged due to collision or drop, please stop using it. Do not reform, disassemble or repair the device by yourself. Before the battery is removed, the device must be powered off. After being removed, the battery should be disposed in a safe manner. This device should be kept away from strong electromagnetic fields (such as induction cookers) when charging. Children are not allowed to play with, clean and maintain the device. Before connecting to the power supply, please check if the voltage indicated on the device matches the local power supply voltage.

The charging cable should not be used outdoors under the sun. Do not use the device when you are taking a bath or shower. Do not share the brush head with others.

Do not place metal objects next to the device while the power is on.

Warning: With small parts inside, the device is not suitable for children under 3 years old. All packaging materials, such as cartons, inner plastic boxes, stickers, outer film, etc. are not part of the device. Children shall not play with them. For the safety of your child, please dispose of the packaging materials mentioned above properly. Warning: The device contains small parts such as the protective cover of the brush head. Avoid them being swallowed by children. When using for the first time or replacing the brush head with a new one, please dispose of the protective cover of the brush head.

Do not use the device beyond the recommended operating temperature of 5°C - 40°C. Do not rinse the toothbrush or any part of the brush head with water above 60°C.

Caution is required for users with sensitive mouth and the following conditions

If you have recently undergone periodontal surgery or are undergoing treatments, or your dentures are not secure, or your gums continue to bleed profusely, or your mouth is sensitive and fragile, please consult a professional dentist before using the device. The device is not suitable for the disabled, mentally ill or those who lack relevant experience and knowledge. Unless used under the supervision of their guardians.

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Use Instructions

Touch the on/off button and the toothbrush starts working.


Long press can switch the brushing modes. There are two modes: Soft and Standard.

When the power runs low, the red light will flash.


Low Battery Prompt: If the red indicator flashes, it means that it needs to be charged to ensure normal use. The red light will flash when the power is below 20%. **How to Charge:** Please gently place the device on the charging base, which has been connected to the power supply. When charging, the white indicator turns on, and the indicator will be constant on after the charging is finished, which usually takes 3 hours. **Important Note:** It is recommended that the charger is energized only when needed. When not needed, please remove the charger. Leaving the toothbrush on the powered charging base for a prolonged period of time may lead to premature depletion of the built-in battery. **Installation of replacement brush head:** Due to abrasion, the brush head needs to be replaced every 2-3 months for hygiene. Please pull out the brush head or insert it into the oscillating shaft by force in the "straight up" and "straight down" directions, and do not rotate it left or right. After installation, there will be a gap between the brush head and the body for the brush head to swing. **Important Note for the Installation of Brush Head:** Since huge force is required to replace the brush head, it is recommended that the replacement should be completed by parents to prevent potential injury caused by improper handling or force by children.

Daily Cleaning and Maintenance


Cleaning the Brush Heads

After each use, please wash the head and brush, and dry it naturally. Do not pull the brush, or it will be deformed and fall off.


Cleaning the Body

Wash the surface of the handle with clean water, wipe it with soft cloth and dry it naturally. Please do not press the rubber seals on the metal axle with sharp objects, or it may cause damages.

Cleaning the Charger

Scrub the stain on the charger with a soft cloth. Before cleaning the charger, make sure to pull out the power supply plug from the socket. Wipe the charger surface with wet cloth. Do not use it before it dries.

Replacement of Brushing Head

It is suggested by professional dentists that the head should be changed every 2-3 months to maintain optimal brushing performance.

IC WARNING

**This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1.This device may not cause interference.

2.This device must accept any interference, including interference that may cause undesired operation of the device.

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
Recommended Tooth Brushing Method

Brushing Time: Dental professionals (dentists) recommend that it takes up to 2 minutes with a traditional manual toothbrush, so we have designed the brushing time for 120 seconds. Brushing Force: Please gently press the toothbrush head against the tooth surface, with no excessive force needed.

FCC WARNING

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

Specifications

Product Name	Oclean Kids Toothbrush for Kids	Battery Capacity	800 mAh (**)	Package Contents	Toothbrush handle * 1, brush head for children * 2, charging base * 1, manual * 1
Type	Oclean Air2	Power Input	5V = 0.5A	Battery Life	30 Days (*)
Maximum Torque Force	180gf.cm	Charging Time	About 3 Hours	Standard	GB 39669-2020
Rated Power	1.8W	Charging Method	Inductive Charging Base		

We are constantly improving our devices. Please note that if the contents in the manual differ slightly from the actual functions, specifications or design of the devices, no notice will be given. Any claims based on the data, diagrams or text descriptions in this manual will not be accepted. For the guidance of user, the pictures in this manual are for reference only. Please refer to the actual devices. The manual contains important information. Please read it carefully before use and keep it.